

TECHNICAL SPECIFICATIONS

6-FAM DNA - Cartridge

Product information

	6-FAM DNA - Cartridge
Volume	100nt
Nucleotides	Deoxyribonucleotides
Bases Quantity	75 x A, C, G, T
5'-Modification	6-FAM dye
Usage	Un-/single-/dual-labeled DNA synthesis
Shelf Life	1 year
Shelf Life after Activation	1 week
Short-Term Storage	2 months at 2-8°C
Long-Term Storage	-20°C
Combination	BHQ-1 Quencher - Chip; 6-FAM Label - Chip; Standard - Chip
Scope of Delivery	1 reagent cartridge
Required	Kilobaser one-XT