

TECHNICAL SPECIFICATIONS

Standard DNA - Cartridge

Product information

	Standard DNA - Cartridge
Volume	150nt
Nucleotides	Deoxyribonucleotides
Bases Quantity	75 x A, C, G, T
5'-Modification	No
Usage	DNA oligo synthesis
Shelf Life	1 year
Shelf Life after Activation	2 weeks
Short-Term Storage	2 months at room temperature
Long-Term Storage	2-8°C
Combination	Standard - Chip
Scope of Delivery	1 reagent cartridge
Required	Kilobaser one or Kilobaser one-XT