

# TECHNICAL SPECIFICATIONS

## OliPure HIC Purification Kit

### Product information

	OliPure HIC Purification Kit
<b>Properties</b>	HIC column based oligonucleotide purification
<b>Throughput</b>	20µL – 120µL sample volume
<b>Min Oligo Length</b>	12nt
<b>Max Oligo Concentration</b>	250 ng/µL
<b>Recovery Rate</b>	75% (±5%)
<b>Shelf Life</b>	1 year
<b>Storage</b>	2-8°C
<b>Application</b>	Purification of oligonucleotides after synthesis or post-synthesis modifications
<b>Scope of Delivery</b>	10 x Spin-down Tubes with HIC column; 1 x Instruction
<b>Required</b>	Microcentrifuge (fixed-angle or swing-out bucket rotors)