

# TECHNICAL SPECIFICATIONS

## Amino Modifier C6 DNA - Cartridge

### Product information

	Amino Modifier C6 DNA - Cartridge
<b>Volume</b>	100nt
<b>Nucleotides</b>	Deoxyribonucleotides
<b>Bases Quantity</b>	75 x A, C, G, T
<b>5'-Modification</b>	Amino Modifier C6
<b>Usage</b>	Un-/single-labeled DNA synthesis
<b>Shelf Life</b>	1 year
<b>Shelf Life after Activation</b>	1 week
<b>Short-Term Storage</b>	2 months at 2-8°C
<b>Long-Term Storage</b>	-20°C
<b>Combination</b>	Standard - Chip; Amino Modifier C6 Label - Chip
<b>Scope of Delivery</b>	1 reagent cartridge
<b>Required</b>	Kilobaser one-XT