

TECHNICAL SPECIFICATIONS

Hexynyl DNA - Cartridge

Product Information

	Hexynyl DNA - Cartridge
Volume	100nt
Nucleotides	Deoxyribonucleotides
Bases Quantity	75 x A, C, G, T
5'-Modification	Hexynyl
Usage	Click chemistry
Shelf Life	1 year
Shelf Life after Activation	1 week
Short-Term Storage	2 months at 2-8°C
Long-Term Storage	-20°C
Combination	Standard - Chip; Hexynyl Label - Chip
Scope of Delivery	1 reagent cartridge
Required	Kilobaser One-XT