

TECHNICAL SPECIFICATIONS

Hexynyl Modification - Chip

Product Information

	Hexynyl Label - Chip
Properties	Single-use; microcolumn incl.
Chip Type	Microfluidic
End Product	5'-hexynyl modified DNA oligo
Material	Polymer
Shelf Life	1 year
Short-term Storage	Dry; room temperature
Long-term storage	Dry; room temperature
Combination	Hexynyl DNA - Cartridge
Scope of Delivery	1 microfluidic chip
Required	Hexynyl DNA - Cartridge; Kilobaser One-XT