

TECHNICAL SPECIFICATIONS

Kilobaser N₂Go

Product Information

	Kilobaser N ₂ Go
Description	Portable Nitrogen Generator
Output Pressure	max. 6.5 bar
Flow Rate	0.8-2 l/min (adjustable)
Internal Gas Storage	2 l
Pressure Regulation	Automatic
Gas Source	Ambient air
Power Input	24 VDC / 8 A
Compatibility	Kilobaser One; Kilobaser One-XT
Scope of Delivery	1 N2Go unit; power supply (24 VDC); pneumatic gas tubing; brief set-up instruction
Instrument Size	400x270x315 mm
Weight	15kg
Warranty	1 year (included); up to 5 additional years (optional)