

# TECHNICAL SPECIFICATIONS

## OliPure Purification Kit

### Product information

	OliPure Purification Kit
<b>Properties</b>	Magnetic beads based oligo purification
<b>Throughput</b>	10 purifications
<b>Min Oligo Length</b>	12nt
<b>Max Oligo Concentration</b>	1nM
<b>Min Elution Volume</b>	10 $\mu$ L
<b>Recovery Rate</b>	>85%
<b>Shelf Life</b>	6 months
<b>Storage</b>	2-8°C; do not freeze
<b>Application</b>	Low-scale oligonucleotide purification after synthesis or post synthesis modification
<b>Scope of Delivery</b>	10 x 50 $\mu$ L Magnetic Beads Suspension; 2 x 500 $\mu$ L Binding Buffer; 2 x 500 $\mu$ L Washing Buffer; 1x Instruction
<b>Required</b>	OliPure Magnetic Rack