

# TECHNICAL SPECIFICATIONS

## 2'-MOE RNA - Cartridge

### Product information and comparison

	2'-MOE RNA - Cartridge
<b>Volume</b>	150nt
<b>Nucleotides</b>	2'-Methoxyethyl Ribonucleotides
<b>Bases Quantity</b>	75 x A, C, G, U
<b>5'-Modification</b>	No
<b>Usage</b>	For modified RNA Synthesis
<b>Shelf Life</b>	1 year
<b>Shelf Life after Activation</b>	2 weeks
<b>Short-Term Storage</b>	2 month at 2-8°C
<b>Long-Term Storage</b>	-20°C
<b>Combination</b>	Standard Chip;
<b>Scope of Delivery</b>	One Reagent Cartridge
<b>Required</b>	Kilobaser one-XT
<b>Price</b>	EUR 750,00 USD 870.00