

TECHNICAL SPECIFICATIONS

2'-OMe RNA - Cartridge

Product information and comparison

	2'-OMe RNA - Cartridge
Volume	150nt
Nucleotides	2'-O-methyl Ribonucleotides
Bases Quantity	75 x A, C, G, U
5'-Modification	No
Usage	For modified RNA Synthesis
Shelf Life	1 year
Shelf Life after Activation	2 weeks
Short-Term Storage	2 month at 2-8°C
Long-Term Storage	-20°C
Combination	Standard Chip;
Scope of Delivery	One Reagent Cartridge
Required	Kilobaser one-XT
Price	EUR 750,00 USD 870.00