

# TECHNICAL SPECIFICATIONS

## 6-FAM DNA - Cartridge

### Product information and comparison

	6-FAM DNA - Cartridge
<b>Volume</b>	100nt
<b>Nucleotides</b>	Deoxyribonucleotides
<b>Bases Quantity</b>	75 x A, C, G, T
<b>5'-Modification</b>	6-FAM Dye
<b>Usage</b>	For not/single/dual labeled DNA oligo Synthesis
<b>Shelf Life</b>	1 year
<b>Shelf Life after Activation</b>	1 week
<b>Short-Term Storage</b>	2 month at 2-8°C
<b>Long-Term Storage</b>	-20°C
<b>Combination</b>	BHQ-1 Chip; 6-FAM Label Chip
<b>Scope of Delivery</b>	One Reagent Cartridge
<b>Required</b>	Kilobaser one-XT
<b>Price</b>	EUR 75,00 USD 100.00