

TECHNICAL SPECIFICATIONS

Amino Modifier C6 DNA - Cartridge

Product information and comparison

	Amino Modifier C6 DNA - Cartridge
Volume	100nt
Nucleotides	Deoxyribonucleotides
Bases Quantity	75 x A, C, G, T
5'-Modification	6-(4-Monomethoxytritylamino)hexyl-(2-cyanoethyl)-(N,N-diisopropyl)-phosphoramidite
Usage	For not/single labeled DNA Synthesis
Shelf Life	1 year
Shelf Life after Activation	1 week
Short-Term Storage	2 month at 2-8°C
Long-Term Storage	-20°C
Combination	Standard Chip; Amino Modifier C6 Label Chip
Scope of Delivery	One Reagent Cartridge
Required	Kilobaser one-XT
Price	EUR 75,00 USD 100.00