

TECHNICAL SPECIFICATIONS

Standard DNA - Cartridge

Product information and comparison

	Standard DNA - Cartridge
Volume	150nt
Nucleotides	Deoxyribonucleotides
Bases Quantity	75 x A, C, G, T
5'-Modification	No
Usage	For DNA Oligo Synthesis
Shelf Life	1 year
Shelf Life after Activation	2 weeks
Short-Term Storage	2 month at room temperature
Long-Term Storage	2-8°C
Combination	Standard Chip;
Scope of Delivery	One Reagent Cartridge
Required	Kilobaser one or Kilobaser one-XT
Price	EUR 75,00 USD 100.00