

TECHNICAL SPECIFICATIONS

6-FAM Label - Chip

Product information and comparison

	6-FAM Label - Chip
Properties	Single-use; Microcolumn incl.
Chip Type	Microfluidic
Usage	For single 6-FAM labeled DNA Oligo Synthesis
Material	Polymer
Shelf life	1 year
Short-term Storage	dry; room temperature
Long-term storage	dry; room temperature
Combination	Kilobaser 6-FAM DNA cartridge
Scope of Delivery	One Microfluidic Chip
Required	Kilobaser one-XT
Price	EUR 65,00 USD 75.00