

TECHNICAL SPECIFICATIONS

Amino Mod C6 Label - Chip

Product information and comparison

	Amino Mod C6 Label - Chip
Properties	Single-use; Microcolumn incl.
Chip Type	Microfluidic
Usage	For 5'-Amino Modification C6 labeled DNA Oligo Synthesis
Material	Polymer
Shelf life	1 year
Short-term Storage	dry; room temperature
Long-term storage	dry; room temperature
Combination	Kilobaser Amino Mod C6 DNA Cartridge
Scope of Delivery	One Microfluidic Chip
Required	Kilobaser one-XT
Price	EUR 30,00 USD 35.00