

TECHNICAL SPECIFICATIONS

Phosphorylation - Chip

Product information and comparison

	Phosphorylation - Chip
Properties	Single-use; Microcolumn incl.
Chip Type	Microfluidic
Usage	For 5'-Phosphorylation labeled DNA Oligo Synthesis
Material	Polymer
Shelf life	1 year
Short-term Storage	dry; room temperature
Long-term storage	dry; room temperature
Combination	Kilobaser Phosphorylation DNA Cartridge
Scope of Delivery	One Microfluidic Chip
Required	Kilobaser one-XT
Price	EUR 20,00 USD 25.00