

TECHNICAL SPECIFICATIONS

Standard - Chip

Product information and comparison

	Standard - Chip
Properties	Single-use; Microcolumn incl.
Chip Type	Microfluidic
Usage	For DNA / RNA Oligo Synthesis
Material	Polymer
Shelf life	1 year
Short-term Storage	dry; room temperature
Long-term storage	dry; room temperature
Combination	All Kilobaser cartridges
Scope of Delivery	One Microfluidic Chip
Required	Kilobaser one or Kilobaser one-XT
Price	EUR 12,50 USD 15.00