

TECHNICAL SPECIFICATIONS

OliPure Kit

Product information and comparison

	OliPure Kit
Properties	Magnetic beads based oligo purification
Throughput	10 Preps
Min Oligo Length	12nt
Max Oligo Concentration	1nM
Min Elution Volume	10 μ L
Recovery Rate	>85%
Shelf life	6 Month
Storage	2-8°C; Do not freeze
Application	Low-scale oligonucleotide purification after synthesis or post synthesis modification
Scope of Delivery	10 x 50 μ L Magnetic Beads Suspension; 2 x 500 μ L Binding Buffer; 2 x 500 μ L Washing Buffer; 1x Instruction
Required	Magnetic Separation Rack
Price	EUR 75,00 USD 100.00